

This document is provided for product evaluation purposes only. It is not intended to be used in place of the package insert shipped with the product.



## DINUCLEOTIDE STANDARD Formamido Cytidine Dimer d(P<sup>f</sup>pC)


This product was manufactured in a facility which has a Quality Management System that is certified as being in compliance with ISO 13485.

ZMC Catalog #: 0801211

The development of this product was supported in part by grants from the National Cancer Institute (R42 CA80451) and the New York Center for Biotechnology.

FOR RESEARCH USE ONLY. Not for *in vitro* Diagnostic Use

### Formamido Cytidine Dimer d(P<sup>f</sup>pC)



d(P<sup>f</sup>pC): Dinucleotide monophosphate where P<sup>f</sup> is the nucleoside bearing the 1'-formamido-2'-deoxy-β-D-ribose modification from a pyrimidine nucleoside.

#### PRODUCT CHARACTERISTICS

d(P<sup>f</sup>pC) is useful as an internal standard for end-labeling as well as LC-MS analysis in the investigation of the nature and consequences of DNA damage. d(P<sup>f</sup>pC) has been chemically synthesized to improve reliability and accuracy and purified by HPLC.<sup>(1)</sup> The modified monophosphate is derived by oxidation of pyrimidines. Identification of the product and an index of its purity was accomplished by proton NMR spectroscopy at 400 MHz.<sup>(2)</sup>


#### RECOMMENDED USAGE

d(P<sup>f</sup>pC) may be used as an internal standard for various quantitations. This product is supplied lyophilized at an OD of 1.0. Reconstitute in water to 1.0 ml.

#### STORAGE AND STABILITY

Lyophilized product is stable for 1 year when stored at 2-8°C. Recommend aliquoting and storing at -20°C when resuspended in water.

## NMR Spectrum



### REFERENCES

- (1) A.E. Maccubbin, H.B. Patrzyc, N. Ersing, Budzinski, J.B. Dawidzik, J.C. Wallace, H. Iijima, H. C. Box, *Biochimica et Biophysica Acta* 1454 (1999) 80-88.
- (2) H.C. Box, H.G. Freund, E.E. Budzinski, J.C. Wallace and A.E. Maccubbin, *Radiation Research*, 141. 91-94 (1995)

ZeptoMetrix Corporation  
872 Main Street  
Buffalo, New York 14202

Office Phone: 716-882-0920  
Fax: 716-882-0959  
Customer Service: 800-274-5487  
[www.zeptometrix.com](http://www.zeptometrix.com)  
e-mail: [support@zeptometrix.com](mailto:support@zeptometrix.com)  
[custserv@zeptometrix.com](mailto:custserv@zeptometrix.com)