



EIA KITS

Enzyme-linked immunosorbent assay (ELISA) is a biochemistry assay method that uses a solid-phase enzyme immunoassay (EIA) to detect the presence of a substance, usually an antigen, in liquid samples.

ELISA tests are commonly carried out using a micro well format and usually use antibodies in the sample linked to an enzyme, which when added to the enzyme's substrate, creates a detectable signal. The detectable signal is usually visualised by a colour change to the substrate.

We provide a range of ELISA kits which are ready to use, with break apart wells. We offer the following range of ELISAs:

- The RETROTEK™ range is designed for the detection and quantitation of retroviral antigens from the retroviruses HIV-1, SIV and HTLV found in cell culture, serum, plasma or other biological fluids
- IMMUNOTEK™ kits can detect and quantify various immunoglobulins from many species including Human, Chicken, Cow, Goat, Rabbit, Rat and Mouse.

We also provide HHV-6 IgG antibody and KSHV/HHV8 IgG antibody ELISA kits.

Product code	Name of the product	Sub-category 2	Product size	Shipping conditions	Webpage
0801060	10x Plate Wash Buffer	Accessory Kits-ELISA Reagents	125 ml	+2° to +8° C.	https://www.h-h-c.com/product/10x-plate-wash-buffer/
0801179	ZeptoCoat	Accessory Kits-ELISA Reagents	60 ml	+2° to +8° C.	https://www.h-h-c.com/product/zeptocoat/
0801184	ZeptoBind	Accessory Kits-ELISA Reagents	60 ml	+2° to +8° C.	https://www.h-h-c.com/product/zeptobind/
0801185	ZeptoBlock	Accessory Kits-ELISA Reagents	120 ml	+2° to +8° C.	https://www.h-h-c.com/product/zeptoblock/
0801188	ELISA Construction System	Accessory Kits-ELISA Reagents	5 plates	+2° to +8° C.	https://www.h-h-c.com/product/elisa-construction-system/
0801179L	ZeptoCoat	Accessory Kits-ELISA Reagents	1 L	+2° to +8° C.	https://www.h-h-c.com/product/zeptocoat-2/
0801184L	ZeptoBind	Accessory Kits-ELISA Reagents	1 L	+2° to +8° C.	https://www.h-h-c.com/product/zeptobind-2/
0801185L	ZeptoBlock	Accessory Kits-ELISA Reagents	1 L	+2° to +8° C.	https://www.h-h-c.com/product/zeptoblock-2/
10-601-001	Immunoassay Stabilizer	Accessory Kits-ELISA Reagents	1 L	+2° to +8° C.	https://www.h-h-c.com/product/immunoassay-stabilizer/
15-401-000	HHV-6 IgG Antibody ELISA Kit	Antibody	Kit	+2° to +8° C.	https://www.h-h-c.com/product/hhv-6-igg-antibody-elisa-kit/
15-501-000	KSHV/HHV-8 IgG Antibody ELISA Kit	Antibody	Kit	+2° to +8° C.	https://www.h-h-c.com/product/kshv/hhv-8-igg-antibody-elisa-kit/
0801009	Rat IgG ELISA	Antigen	96 det.	+2° to +8° C.	https://www.h-h-c.com/product/rat-igg-elisa/
0801010	Chicken IgY ELISA	Antigen	96 det.	+2° to +8° C.	https://www.h-h-c.com/product/chicken-igy-elisa/
0801180	Mouse IgG ELISA	Antigen	96 det.	+2° to +8° C.	https://www.h-h-c.com/product/mouse-igg-elisa/

Product code	Name of the product	Sub-category 2	Product size	Shipping conditions	Webpage
0801181	Mouse IgM ELISA	Antigen	96 det.	+2° to +8° C.	https://www.h-h-c.com/product/mouse-igm-elisa/
0801182	Human IgG ELISA	Antigen	96 det.	+2° to +8° C.	https://www.h-h-c.com/product/human-igg-elisa/
0801183	Human IgM ELISA	Antigen	96 det.	+2° to +8° C.	https://www.h-h-c.com/product/human-igm-elisa/
0801197	Human IgA ELISA	Antigen	96 det.	+2° to +8° C.	https://www.h-h-c.com/product/human-iga-elisa/
0801198	Bovine IgG ELISA	Antigen	96 det.	+2° to +8° C.	https://www.h-h-c.com/product/bovine-igg-elisa/
0801202	Rabbit IgG ELISA	Antigen	96 det.	+2° to +8° C.	https://www.h-h-c.com/product/rabbit-igg-elisa/
0801221	Horse IgG ELISA	Antigen	96 det.	+2° to +8° C.	https://www.h-h-c.com/product/horse-igg-elisa/
0801223	Goat IgG ELISA	Antigen	96 det.	+2° to +8° C.	https://www.h-h-c.com/product/goat-igg-elisa/
0801002	HIV-1 p24 Antigen ELISA 2.0 Kit	Antigen	96 det.	+2° to +8° C.	https://www.h-h-c.com/product/hiv-1-p24-antigen-elisa-2.0-kit/
0801008	HIV-1 p24 Antigen ELISA 2.0 Kit	Antigen	480 det.	+2° to +8° C.	https://www.h-h-c.com/product/hiv-1-p24-antigen-elisa-2.0-kit-2/
0801096	HIV-1 p24 ICx - CRx Kit	Antigen	100 det.	+2° to +8° C.	https://www.h-h-c.com/product/hiv-1-p24-icx---crx-kit/
0801111	HIV-1 p24 Antigen ELISA Kit	Antigen	96 det.	+2° to +8° C.	https://www.h-h-c.com/product/hiv-1-p24-antigen-elisa-kit/
0801116	HTLV p19 Antigen ELISA Kit	Antigen	96 det.	+2° to +8° C.	https://www.h-h-c.com/product/htlv-p19-antigen-elisa-kit/
0801137	HIV-1 p24 Extended Range Kit	Antigen	kit	+2° to +8° C.	https://www.h-h-c.com/product/hiv-1-p24-extended-range-kit/
0801169	SIV p27 Antigen ELISA Kit	Antigen	96 det.	+2° to +8° C.	https://www.h-h-c.com/product/siv-p27-antigen-elisa-kit/
0801200	HIV-1 p24 Antigen ELISA Kit	Antigen	480 det.	+2° to +8° C.	https://www.h-h-c.com/product/hiv-1-p24-antigen-elisa-kit-2/
0801201	SIV p27 Antigen ELISA Kit	Antigen	480 det.	+2° to +8° C.	https://www.h-h-c.com/product/siv-p27-antigen-elisa-kit-2/

Helvetica Health Care supplies Life Science products to manage health and improve the quality of life. We are committed to unrivalled service with 48 hour, risk free delivery and order tracking from despatch to receipt.

Visit us and get in touch



www.h-h-c.com