

HIV DATA SUMMARY

PANEL NUMBER	DONOR NUMBER	PATIENT DRAWS	DAYS SPANNED	HIV RNA	HIV-1 bDNA	HIV Ag	ANTI-HIV	HIV AB 1/2+O	HIV-Ag/Ab COMBO
HIV Pre/Post									
HIV10234	67101	3	39		*	POS	NEG	*	*
Anti-HIV-1 False Positive									
HIV9027	65633	9	29	*	POS	NEG	POS	POS	POS
HIV9078	79065	5	70	*	*	POS	NEG	NEG	NEG
Anti-HIV-1									
HIV6243	62238	10	35	POS	*	POS	POS	POS	POS
HIV6247	63602	10	31	POS	POS	POS	POS	POS	POS
HIV6248	63331	7	26	POS	POS	POS	POS	POS	POS
HIV9011	64954	11	39	*	POS	POS	POS	POS	POS
HIV9012	65389	8	23	*	POS	POS	POS	POS	POS
HIV9013	65404	7	26	*	POS	POS	NEG	NEG	POS
HIV9014	65522	7	31	*	POS	POS	POS	POS	POS
HIV9015	65661	8	38	*	POS	POS	POS	POS	POS
HIV9016	65790	10	35	*	POS	POS	NEG	NEG	POS
HIV9018	66575	11	36	*	POS	POS	POS	POS	POS
HIV9019	65685	3	38	*	POS	POS	POS	POS	POS
HIV9020	67512	22	98	*	POS	POS	POS	POS	POS
HIV9021	67485	17	58	*	POS	POS	POS	POS	POS
HIV9022	64578	9	33	*	POS	POS	POS	POS	POS
HIV9023	67706	22	86	*	POS	POS	POS	NEG	POS
HIV9024	64469	12	54	*	POS	POS	POS	NEG	POS
HIV9025	67996	12	92	*	POS	POS	POS	POS	POS
HIV9026	68205	7	45	*	POS	POS	POS	POS	POS
HIV9028	67860	7	56	*	POS	POS	NEG	NEG	POS
HIV9029	68262	12	46	*	POS	POS	POS	NEG	POS
HIV9030	68582	16	55	*	POS	POS	POS	POS	POS
HIV9031	66048	19	158	*	POS	POS	POS	POS	POS
HIV9032	68106	14	56	*	POS	POS	POS	POS	POS
HIV9033	66686	16	85	*	POS	POS	POS	POS	POS
HIV9034	66632	13	54	*	POS	POS	POS	POS	POS
HIV9075	62216	5	48	POS	*	POS	POS	POS	POS
HIV9076	63753	9	74	POS	*	POS	POS	POS	POS
HIV9077	75018	27	92	*	*	POS	POS	POS	POS
HIV9079	75062	25	98	POS	*	POS	POS	POS	POS
HIV9081	63215	4	34	POS	*	POS	POS	POS	POS
HIV9082	77840	4	21	POS	*	POS	POS	*	*
HIV9084	77600	4	49	POS	POS	NEG	POS	*	*
HIV9089	65376	6	24	*	POS	POS	POS	POS	POS
HIV9096	1043487	6	22	POS	*	*	POS	*	*
HIV12007	73695	9	134	POS	*	NEG	POS	POS	POS
HIV12008	73698	13	43	POS	*	POS	POS	POS	POS
