

A

Adult Human Elutriated Lymphocytes – Frozen

Catalog Number: 07-212-001

Lot Number: C0618

Product Description: Frozen viable adult human lymphocytes elutriated from source leukocytes.

Unit Size: 1.0 ml

Blood Donor: Race: Caucasian
Sex: Male
Blood Type: Unknown

Freeze Medium: RPMI with 20% fetal bovine serum, gentamicin, and 10% DMSO.

QUALITY CONTROL DATA

Sterility Tests: Bacteria, yeasts, and fungi were not detected (USP-NF).
Mycoplasma was not detected by MycoAlert™ Mycoplasma Detection Kit (Lonza, LT07).

Pre-Freeze Cell Count: 4×10^7 cells/ml

Post-Thaw Cell Count: 3×10^7 cells/ml

PRODUCT DETAILS

Shipping and Storage: This product is shipped frozen on dry ice. **Store at -130°C or below upon receipt.**

Safe Handling Recommendation: This product was tested and found negative for HBsAg, anti-HIV-1/2, anti-HCV, and HIV antigen. No test guarantees a product to be completely non-infectious. Therefore, all human blood derived material should be considered as potentially infectious and handled in accordance with biosafety guidelines defined in the BMBL, NIH-CDC HHS Publication No. 93-8395.

Product Disclaimer: Product variances may result from donor characteristics, including overall donor health, blood draw timing, and other conditions. Variances may include cell counts and cell viability and are **NOT** the responsibility of ABI. Variation in laboratory technique, methods, equipment, and testing materials (i.e., Human AB sera, cytokines, etc.), may create unexpected differences in individual lot of cells' growth performance.

This product is for research use only.
Not for use in diagnostic procedures.



Quality Control

08/30/2012

Date