

MONOBODIES

ANTI-HIV-1 p24

Clone 39/5.4A

ZMC Catalog #: 0801080 **Unit Size:** 100 µg

ZMC Catalog #: 0801136 **Unit Size:** 1 mg

PRODUCT CHARACTERISTICS

This murine monoclonal antibody reacts with Human Immunodeficiency Virus Type 1(HIV-1) p24 protein. The antibody is of the IgG₁ subclass and has been purified from serum-free culture supernatant by protein G chromatography. The antibody is supplied as purified immunoglobulin at a concentration of approximately 1 mg/ml in phosphate-buffered saline. MONOBODIES Anti-HIV-1 p24 (clone 39/5.4A) exhibits reactivity with HIV-1 infected cultures using indirect immunofluorescence. It also exhibits reactivity with viral lysates examined using Western blotting. In addition, the antibody immunoprecipitates p24 protein from radiolabeled infected cell extracts. This antibody may be useful as a solid phase absorbent in an ELISA format.

RECOMMENDED USAGE

Antibody dilutions should be prepared using buffers containing suitable protein in order to stabilize antibody activity. The antibody is used to detect HIV-1 core protein in virus infected cells as well as in viral or infected cell lysates. Studies on core antigen synthesis and metabolism can be performed using Western blotting or radioimmunoprecipitation analysis.

STORAGE

Aliquot and store at -20°C.

1004159

Rev.: 04/00

All products are for Research Use Only and not for *in vitro* Diagnostic Use

This Package Insert is provided for product evaluation purposes only and is not intended to be used in place of the Package Insert shipped with the product.