

***Rhinovirus A16
Culture Fluid (Heat Inactivated) (1 mL)**

Catalog Number: 0810285CFHI

PRODUCT DESCRIPTION:

Rhinovirus A16 consists of a non-enveloped capsid containing a positive-sense single-stranded RNA.

Each frozen aliquot contains 1 mL of heat inactivated viral culture fluid. The pre-inactivation titer was determined from an infectious aliquot.

Viral inactivation is verified after heat inactivation by the absence of viral growth in tissue culture-based infectivity assays.

INTENDED USE:

Heat inactivated viral culture fluids are sold as consumable testing materials. The suitability and performance characteristics should be determined by your laboratory for each intended usage.

These products are NOT intended for use in the manufacture or processing of injectable products subject to licensure under section 351 of the Public Health Service Act or for any other product intended for administration to humans.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

BIOSAFETY:

Rhinovirus is a Biosafety Level 2 microorganism when infectious and should be used within a Biological Safety Level 2 facility or cabinet. Please consult your institution's regulations regarding the use of this product. For a detailed discussion on biological safety see the 5th edition of Biosafety in Microbiological and Biomedical Laboratories (BMBL), published by the CDC at <http://www.cdc.gov/biosafety/publications/bmb15/index.htm>.

PRECAUTIONS:

- Use Universal Precautions, this product is **potentially biohazardous**.
- Repetitive freezing and thawing is not recommended (aliquot material if necessary).
- To avoid cross-contamination, use separate pipette tips for all reagents.









RECOMMENDED STORAGE:

Heat inactivated viral culture fluids should be stored at -65°C or below.

***This strain is sourced and used under license from the National Collection of Pathogenic Viruses (NCPV®), Public Health England.**

®Registered trademarks are the property of their respective owners.

This product was manufactured in a facility which has a Quality Management System that is ISO 13485 certified.

 REF	Catalog Number		Temperature Limitation
 LOT	Lot Number		Expiration Date
 RUO	For Research Use Only		Biological Risk
	Manufacturer		Date of Manufacture

ZeptoMetrix Corporation
878 Main Street
Buffalo, NY 14202
800-274-5487
www.zeptometrix.com