

**SAR SARS-Related Coronavirus 2 (SARS-CoV-2)  
Lineage B.1.617.2; Delta Variant\*  
Isolate: USA/PHC658/2021  
Culture Fluid (UV Inactivated) (1 mL)**

Catalog Number:  
0810624UV

## PRODUCT DESCRIPTION:

SARS-Related Coronavirus 2 Lineage B.1.617.2; Delta Variant (Isolate: USA/PHC658/2021) is an enveloped, positive-sense single-stranded RNA virus from the *Coronaviridae* family and the *Betacoronaviridae* genus. Severe acute respiratory syndrome-related coronavirus 2 (SARS-CoV-2), isolate hCoV-19/USA/PHC658/2021 was isolated in Memphis, Tennessee, USA in May 2021. Under the nomenclature system introduced by GISAID (Global Initiative on Sharing All Influenza Data), SARS-CoV-2, isolate hCoV-19/USA/PHC658/2021 is assigned lineage B.1.617.2 using Phylogenetic Assignment of Named Global Outbreak LINEages (PANGOLIN) tool. SARS-CoV-2, lineage B.1.617.2 was first detected in India and labelled as a variant of concern (VOC). It was labelled as Delta variant by the World Health Organization (WHO).

Each frozen aliquot contains 1 mL of UV inactivated, clarified viral culture fluid. The pre-inactivation titer was determined from an infectious aliquot.

Viral inactivation is verified after UV inactivation by the absence of viral growth in tissue culture-based infectivity assays.

## INTENDED USE:

UV inactivated viral culture fluids are sold as consumable testing materials. The suitability and performance characteristics should be determined by your laboratory for each intended usage.

These products are NOT intended for use in the manufacture or processing of injectable products subject to licensure under section 351 of the Public Health Service Act or for any other product intended for administration to humans.

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

## BIOSAFETY:

Inactivated SARS-CoV-2 should be used within a Biological Safety Level 2 facility or biosafety cabinet. Please consult your institution's regulations regarding the use of this product. For a detailed discussion on biological safety see the 5<sup>th</sup> edition of Biosafety in Microbiological and Biomedical Laboratories (BMBL), published by the CDC at <http://www.cdc.gov/biosafety/publications/bmbl5/index.htm>

## PRECAUTIONS:

- Use Universal Precautions, this product is **potentially biohazardous**.
- Repetitive freezing and thawing is not recommended (aliquot material if necessary).
- To avoid cross-contamination, use separate pipette tips for all reagents.

## RECOMMENDED STORAGE:

UV inactivated viral culture fluids should be stored at -65°C or below.

**\*The following reagent was obtained through BEI Resources, NIAID, NIH: SARS-Related Coronavirus 2, Isolate hCoV-19/USA/PHC658/2021 (Lineage B.1.617.2; Delta Variant), NR-55611, contributed by Dr. Richard Webby and Dr. Anami Patel.**

PI0810624UV  
Revision: 00  
Effective Date: 08/02/2021

	Catalog Number		Temperature Limitation
	Batch Code		Expiration Date
	For Research Use Only		Biological Risk
	Manufacturer		

[www.ZeptoMetrix.com](http://www.ZeptoMetrix.com)

ZeptoMetrix LLC • 878 Main Street, Buffalo, NY 14202 USA • Tel (800) 274-5487

This product was manufactured in a facility which has a Quality Management System that is ISO 13485 certified.

PCA# N/A  
Page 1 of 1