

**A****KSHV/HHV-8 IgG Antibody ELISA Kit**

Catalog Number: 15-501-000

Lot Number: F0202

Product Description: Anti-Kaposi's Sarcoma-Associated Herpesvirus/Human Herpesvirus 8 Immunoglobulin G Enzyme-Linked Immunosorbent Assay Kit.

Expiration Date: 04-2017

Components	Catalog Number	Lot Number
A. KSHV/HHV-8 Antigen Coated ELISA Plate (12 x 8 well strips)	15-503-000	F0202-PL
B. Specimen Diluent (4X)	15-507-000	F0202-SD
C. Negative Control	15-509-000	F0202-N
D. Positive Control	15-508-000	F0202-P
E. ELISA Wash Buffer (20X)	15-208-000	F0202-W
F. HRP Conjugate	15-506-000	F0202-C
G. Conjugate Substrate	15-204-000	F0202-CS
H. Stop Solution	15-207-000	F0202-S
I. ELISA Accessory Pack	15-504-000	F0202-A
J. Package Insert (1 copy)	15-501-000	N/A (revised 04/2015)

**QUALITY CONTROL DATA**

Manufacturing Parameters:	Specifications:	Results:
<b>Plates:</b> Packing:	Sealed pouch and desiccant	Pass
<b>Specimen Diluent (4X):</b> Visual Inspection: pH (1:4 dilution):	Clear Liquid 7.1-7.6	Pass 7.4
<b>Negative Control:</b> Reactivity: Visual Inspection:	OD 450-650 (nm) ≤ 0.100 Clear Liquid	0.021 Pass
<b>Positive Control:</b> Reactivity: Visual Inspection:	OD 450-650 (nm) ≥ 0.300 Clear Liquid	2.185 Pass
<b>ELISA Wash Buffer (20X):</b> Precipitation: pH (1:20 dilution): Osmolality (1:20 dilution):	None 7.20-7.40 250-350 mOsm/kg	Pass 7.4 300
<b>HRP Conjugate:</b> Visual Inspection:	Clear Liquid	Pass
<b>Conjugate Substrate:</b> Visual Inspection:	Clear Liquid	Pass
<b>Stop Solution:</b> Visual Inspection:	Clear Liquid	Pass
<b>Kit Performance:</b>	Pass	Pass

**PRODUCT DETAILS**

Shipping and Storage: This product is shipped at 2-8°C.  
Store at 2-8°C upon receipt.

This product is for research use only.  
Not for use in diagnostic procedures.

Quality Control

04-24-2015

Date

Quick. Easy. **Technical support.**  
[www.abionline.com/contact-us/](http://www.abionline.com/contact-us/)



**CERTIFICATE OF ANALYSIS**