



HCV Seroconversion Panel Donor No. 22708D
Part Number: HCV9094

Donor Demographics: Donor source plasma was collected from a U.S. licensed blood collection facility.
Donor was un-treated and asymptomatic throughout the collection period.

| PANEL MEMBER | DATE OF DRAW | HCV | QUEST HCV Ab | ROCHE HCV Ampliscreen | ROCHE HCV Monitor | ORTHO ECi HBsAg | ORTHO ECi ANTI-HCV | ORTHO ECi ANTI-HIV ½ |
|--------------|--------------|----------------|----------------------|-----------------------|----------------------|-----------------------|-----------------------|-----------------------|
| | | TEST DATE: UNK | TEST DATE: 2/25/2010 | TEST DATE: 4/19/2010 | TEST DATE: 4/20/2010 | TEST DATE: 11/16/2018 | TEST DATE: 11/16/2018 | TEST DATE: 11/16/2018 |
| 9094-01 | 05/03/2004 | NEGATIVE | 0.28 | NEGATIVE | NOT DETECTED | 0.06 | 0.22 | 0.14 |
| 9094-02 | 05/05/2004 | NEGATIVE | 0.34 | NEGATIVE | NOT DETECTED | 0.06 | 0.18 | 0.15 |
| 9094-03 | 05/10/2004 | NEGATIVE | 1.86 | NEGATIVE | NOT DETECTED | 0.06 | 1.14 | 0.15 |
| 9094-04 | 05/12/2004 | POSITIVE | 2.32 | NEGATIVE | NOT DETECTED | 0.05 | 2.31 | 0.16 |
| 9094-05 | 05/17/2004 | POSITIVE | 3.00 | NEGATIVE | NOT DETECTED | 0.05 | 2.00 | 0.16 |

Notes on Data:

1. This panel is for **Research Use Only. Not for use in diagnostic procedures.**
2. Data is generated at ZeptoMetrix™ LLC and by our Laboratory Partners in the United States and Europe.
3. All specimens collected in this longitudinal series are unadulterated 4% sodium citrated plasma samples collected from a single donor in the United States.
4. No preservatives have been added.
5. All material has been tested and found negative for HBsAg, Anti-HIV1/2, HIV p24 Ag by FDA cleared tests.

ZeptoMetrix LLC • 25 Kenwood Circle, Franklin, MA 02038 • Tel (508) 553-5800 • Fax (508) 520-1525
This Product was manufactured in a facility which has a Quality Management System that is ISO 13485 certified.
For Customer Support, please visit www.zeptometrix.com or email zepto.customerservice@antylia.com

PIHCV9094 Rev02
Effective Date: 10/11/2021