

Baby System®

Portable Graf's cradle

Ultrasound Positioning Device with arm guide for ultrasound probes - Medical speciality : medical sonographer/doctors (orthopedists, pediatricians)/nurses

The Portable Graf's Cradle is a lighter version of Graf's Cradle product code 5003. The new version of the cradle is smaller and lighter and easy to transport but still allows the scan of new born hip with the same efficiency. It includes the positioning device for holding and positioning the newborn and the arm guide for the ultrasound scan probe. The arm guide can be positioned with a magnet (and secured with screws) on the right or left side of the positioning device depending on operator requirements. The portable Graf's Cradle is made with same quality products and attention to newborn safety and operators professional requirements.



DESIGNED, CERTIFIED AND AUTHORIZED BY Dr. Prof. Reinhard Graf



Additional Information

- Product code 5002 - HS 94029000 - size 44x38x22 cms - weight 4.8 KGS
- Alternative product : Graf's cradle product code 5003

[CHECK PRICE ON LINE](https://www.h-h-c.com)