

NATtrol™ Women's Health & STI Products

Molecular Controls & Panels for Infectious Disease Testing

ZeptoMetrix® NATtrol™ Women's Health & STI products are the most trusted brand in third party molecular diagnostic QC.

Utilizing ZeptoMetrix [NATtrol Molecular Controls](#) on a consistent basis helps to monitor testing performance and provides additional confidence in the quality and reliability of clinical laboratory operations.

Refrigerated ready-to-use | Safe, non-infectious | Non-synthetic whole microorganisms | Mimics a true patient sample

Part Number	Product Description	Use*	Platform(s)**
BV			
NATBVNEG-BD	Bacterial vaginosis (BV) Negative Control (6 x 0.15 mL)	RUO	BD Max™
NATBVPOS-BD	Bacterial vaginosis (BV) Positive Control (6 x 0.15 mL)	RUO	BD Max™
CT/NG			
NATCT(434)-ERCL	Chlamydia trachomatis LGV II 434, External Run Control, Low (6 x 1 mL) <i>Target concentration: 5,000 copies/mL</i>	RUO	Real-time PCR thermal cyclers
NATCT(434)-ERCM	Chlamydia trachomatis LGV II 434, External Run Control, Medium (6 x 1 mL) <i>Target concentration: 50,000 copies/mL</i>	RUO	Real-time PCR thermal cyclers
NATCT(434)-6MC	Chlamydia trachomatis Positive Control (6 x 1 mL)	IVD	Cepheid® GeneXpert®
NATCT(434)-6MC-IVD	Chlamydia trachomatis Positive Control (6 x 1.25 mL)	CE IVD	Cepheid® GeneXpert®, Hologic® Panther Fusion®
NATCT(D-UW3)-ERCL	Chlamydia trachomatis serotype D, External Run Control, Low (6 x 1 mL) <i>Target concentration: 5,000 copies/mL</i>	RUO	Real-time PCR thermal cyclers
NATCT(D-UW3)-ERCM	Chlamydia trachomatis serotype D, External Run Control, Medium (6 x 1 mL) <i>Target concentration: 50,000 copies/mL</i>	RUO	Real-time PCR thermal cyclers
NATCT/NGNEG-6MC	Chlamydia trachomatis/Neisseria gonorrhoeae Negative Control (6 x 1 mL)	RUO	Cepheid® GeneXpert®
NATCT/NGNEG-6MC-IVD	Chlamydia trachomatis/Neisseria gonorrhoeae Negative Control (6 x 1.25 mL)	CE IVD	Cepheid® GeneXpert®, Hologic® Panther Fusion®
NATCT/NGP-C	Chlamydia trachomatis/Neisseria gonorrhoeae Panel (17 x 1.2 mL)	RUO	Cepheid® GeneXpert®
NATNG-ERCL	Neisseria gonorrhoeae Z017, External Run Control, Low (6 x 1 mL) <i>Target concentration: 5,000 copies/mL</i>	RUO	Real-time PCR thermal cyclers
NATNG-ERCM	Neisseria gonorrhoeae Z017, External Run Control, Medium (6 x 1 mL) <i>Target concentration: 50,000 copies/mL</i>	RUO	Real-time PCR thermal cyclers
NATNG-6MC	Neisseria gonorrhoeae Positive Control (6 x 1 mL)	RUO	Cepheid® GeneXpert®
NATNG-6MC-IVD	Neisseria gonorrhoeae Positive Control (6 x 1.25 mL)	CE IVD	Cepheid® GeneXpert®
CT/NG/TV			
NATCTVPOS-BD	Candida/TV Positive Control (6 x 0.15 mL)	RUO	BD Max™
NATCTNGTV-NEG	CT/NG/TV Negative Control (6 x 1.0 mL)	RUO	Visby Medical™
NATCTNGTV-POS	CT/NG/TV Positive Control (6 x 1.0 mL)	RUO	Visby Medical™
NATTVNEG-6MC	Trichomonas vaginalis Negative Control (6 x 1.2 mL)	RUO	Cepheid® GeneXpert®
NATTVPOS-6MC	Trichomonas vaginalis Positive Control (6 x 1.2 mL)	RUO	Cepheid® GeneXpert®
NATTVGP-C	Trichomonas vaginalis Verification Panel (17 x 0.7 mL)	RUO	Cepheid® GeneXpert®
NATVNEG-6C NEW	Vaginal Negative Control (6 x 0.9 mL)	IVD	Cepheid® GeneXpert®
NATVPOS-6C NEW	Vaginal Positive Control (6 x 0.9 mL)	IVD	Cepheid® GeneXpert®
NATVP-BD	Vaginal Panel (24 x 0.5 mL)	RUO	BD Max™
GBS			
NATLAC-6MC	GBS Negative Control (6 x 0.5 mL)	RUO	BD Max™, Cepheid® GeneXpert®
NATLAC-6MC-IVD	GBS Negative Control (6 x 0.5 mL)	CE IVD	Cepheid® GeneXpert®
NATSAG-6L	GBS Positive Control (6 x 0.25 mL)	RUO	Luminex® ARIES®
NATSAG-6MC	GBS Positive Control (6 x 0.5 mL)	RUO	BD Max™, Cepheid® GeneXpert®
NATSAG-6MC-IVD	GBS Positive Control (6 x 0.5 mL)	CE IVD	Cepheid® GeneXpert®
MGN			
NATMGN-ERC	Mycoplasma genitalium External Run Control (6 x 0.25 mL) <i>Target concentration Ct range: 28–31</i>	RUO	Hologic® Panther Fusion®

*RUO = Research Use Only; IVD = In Vitro Diagnostic; CE IVD = In Vitro Diagnostic with European CE mark

**ZeptoMetrix External Quality Controls and Verification Panels have demonstrated utility on multiple assays. Based upon product configuration, including fill volumes and concentrations of analytes, we can confidently recommend the optimal choice of controls and panels for the indicated platform.

Distributor :
Helvetica Health Care Sàrl
 Voie-de-traz 20
 1218 Grand-Sconnex
 Switzerland

CONTACT US through
our new e-commerce : <https://www.h-h-c.com>
e-mail : incoming@h-h-c.com
Tel. +41 22 5349703

HHC
 Helvetica Health Care