



# NATtrol™

IVD MOLECULAR CONTROLS

## FOR THOSE WHO MUST GET IT RIGHT

ZeptoMetrix is your end-to-end partner for high-quality infectious disease controls for your IVD diagnostic assays. Our ZeptoMetrix NATtrol™ is the most trusted brand of third-party molecular diagnostic QC for infectious disease testing. Using NATtrol Molecular Controls on a consistent basis helps to monitor testing performance and provides additional confidence in the quality and reliability of clinical laboratory operations.

**What really matters to your lab: Refrigerated ready-to-use + Non-infectious + Whole microorganisms + Mimics a true clinical specimen**

- Can be used across various molecular testing platforms and assays
- Available in different formats (single or multiple analytes)
- Full process controls that monitor extraction and amplification steps
- Train and monitor laboratory personnel
- Evaluate lot-to-lot consistency of test kits and assay reagents



**Web**

[www.h-h-c.com](http://www.h-h-c.com)



**Phone**

+41 22 534 97 03



**Email**

[incoming@h-h-c.com](mailto:incoming@h-h-c.com)



**Fax**

+41 22 594 85 07



**Address**

Helvetica Health Care Sàrl – PO Box 1158 - CH – 1211 Geneva 5 AP - Switzerland

## Respiratory

Part Number	Product Description	Use*	Platform(s)**
NATFLUAB-6C-IVD	Influenza A/B Positive Control (6 x 0.5 mL)	CE IVD	Cepheid® GeneXpert®
NATCXVA9-6C-IVD	Influenza/RSV Negative Control (6 x 0.5 mL)	CE IVD	Cepheid® GeneXpert®
NATCV9-6C-IVD	Influenza/RSV/SARS-CoV-2 Negative Control (6 x 0.5 mL)	CE IVD	Cepheid® GeneXpert®, Roche cobas® Liat®
NATFRC-6C-IVD	Influenza/RSV/SARS-CoV-2 Positive Control (6 x 0.5 mL)	CE IVD	Cepheid® GeneXpert®, Roche cobas® Liat®
NATRSV-6C-IVD	Respiratory Syncytial Virus Positive Control (6 x 0.5 mL)	CE IVD	Cepheid® GeneXpert®
NATSARS(COV2)-NEG-IVD	SARS-CoV-2 Negative Control (6 x 0.5 mL)	CE IVD	Abbott ID NOW™, BioFire® FilmArray®, Cepheid® GeneXpert®, GenMark® ePlex®, Luminex® ARIES®, QIAGEN® QIAstat-Dx®, Real-time PCR
NATSARS(COV2)-NEG1-IVD	SARS-CoV-2 Negative Control (6 x 1 mL)	CE IVD	Abbott ID NOW™, BD MAX™, BioFire® FilmArray®, Cepheid® GeneXpert®, GenMark® ePlex®, Luminex® ARIES®, QIAGEN® QIAstat-Dx®, Real-time PCR, Visby Medical™
NATSARS(COV2)-ERC-IVD	SARS-CoV-2 Positive Control (6 x 0.5 mL)	CE IVD	Abbott ID NOW™, BioFire® FilmArray®, Cepheid® GeneXpert®, GenMark® ePlex®, Luminex® ARIES®, QIAGEN® QIAstat-Dx®, Real-time PCR
NATSARS(COV2)-ERC1-IVD	SARS-CoV-2 Positive Control (6 x 1 mL)	CE IVD	Abbott ID NOW™, BD MAX™, BioFire® FilmArray®, Cepheid® GeneXpert®, GenMark® ePlex®, Luminex® ARIES®, QIAGEN® QIAstat-Dx®, Real-time PCR, Visby Medical™

## Gastrointestinal (GI)

Part Number	Product Description	Use*	Platform(s)**
NATCSO-6MC-IVD	C. difficile Negative Control (6 x 0.5 mL)	CE IVD	Cepheid® GeneXpert®
NATCDI-6MC-IVD	C. difficile Positive Control (6 x 0.5 mL)	CE IVD	Cepheid® GeneXpert®
NATNOV-6MC-IVD	Norovirus GI/GII Positive Control 6 X 0.5 mL	CE IVD	Cepheid® GeneXpert®
NATROTA-6MC-IVD	Norovirus Negative Control 6 X 0.5 mL	CE IVD	Cepheid® GeneXpert®

## Women's Health & STI

Part Number	Product Description	Use*	Platform(s)**
NATCT(434)-6MC	C. trachomatis Positive Control (6 x 1 mL)	IVD	Cepheid® GeneXpert®
NATCT(434)-6MC-IVD	C. trachomatis Positive Control Pack (6 x 1.25 mL)	CE IVD	Cepheid® GeneXpert®, Hologic® Panther Fusion®
NATCTNGNEG-6MC	C. trachomatis/N. gonorrhoeae Negative Control (6 x 1 mL)	IVD	Cepheid® GeneXpert®
NATCT-NGNEG-6MC-IVD	C. trachomatis/N. gonorrhoeae Negative Control (6 x 1.25 mL)	CE IVD	Cepheid® GeneXpert®, Hologic® Panther Fusion®
NATCTNGTV-NEG-IVD	CT/NG/TV Negative Control (6 x 1.0 mL)	IVD	Visby Medical™
NATCTNGTV-POS-IVD	CT/NG/TV Positive Control (6 x 1.0 mL)	IVD	Visby Medical™
NATLAC-6MC-IVD	GBS Negative Control (6 x 0.5 mL)	CE IVD	Cepheid® GeneXpert®
NATSAG-6MC-IVD	GBS Positive Control (6 x 0.5 mL)	CE IVD	Cepheid® GeneXpert®
NATNG-6MC	N. gonorrhoeae Positive Control (6 x 1 mL)	IVD	Cepheid® GeneXpert®
NATNG-6MC-IVD	N. gonorrhoeae Positive Control (6 x 1.25 mL)	CE IVD	Cepheid® GeneXpert®
NATTVNEG-6MC-IVD	T. vaginalis Negative Control (6 x 1.2 mL)	CE IVD	Cepheid® GeneXpert®
NATTVPOS-6MC-IVD	T. vaginalis Positive Control (6 x 1.2 mL)	CE IVD	Cepheid® GeneXpert®
NATVNEG-6C	Vaginal Negative Control (6 x 0.9 mL)	IVD	Cepheid® GeneXpert®
NATVPOS-6C	Vaginal Positive Control (6 x 0.9 mL)	IVD	Cepheid® GeneXpert®

## Healthcare-Associated Infections

Part Number	Product Description	Use*	Platform(s)**
NATCSO-6MC-IVD	C. difficile Negative Control (6 x 0.5 mL)	CE IVD	Cepheid® GeneXpert®
NATCDI-6MC-IVD	C. difficile Positive Control (6 x 0.5 mL)	CE IVD	Cepheid® GeneXpert®
NATMRS-6MC-IVD	MRSA Positive Control (6 x 0.5 mL)	CE IVD	Cepheid® GeneXpert®
NATMSSE-6MC-IVD	MRSA/SA Negative Control (6 x 0.5 mL)	CE IVD	Cepheid® GeneXpert®
NATMSSA-6MC-IVD	S. aureus (MSSA) Positive Control (6 x 0.5 mL)	CE IVD	Cepheid® GeneXpert®

\*IVD = In Vitro Diagnostic; CE IVD = In Vitro Diagnostic with European CE mark, Notified Body code 2797

**Web**  
www.h-h-c.com

**Phone**  
+41 22 534 97 03

**Email**  
incoming@h-h-c.com

**Fax**  
+41 22 594 85 07

**Address**  
Helvetica Health Care Sàrl – PO Box 1158 – CH – 1211 Geneva 5 AP – Switzerland